

Rev. 0 - 2020

### KEY FEATURES

For the daily cleaning and disinfection of all types of aspirator/suction systems. Can be used on wet line, dry line and semi wet/dry line systems

- Non foaming & non abrasive
- Contact time from only 15 minutes
- Low hazard classification
- Alkaline based
- Wide effective spectrum
- Blue in colour
- Concentrate solution
- Economical 2% dilution ratio

Viroxid® Aspirator is free from chlorides, aldehyde and phenols. The unique non foaming formula ensures equipment remains in top condition. The highly effective solution is Bactericidal, Yeasticidal, Fungicidal & Tuberculocidal. Also kills Hep B, Hep C and HIV viruses. A concentrated disinfectant solution with an economical dilution ratio of only 2% ensures excellent value per litre of usable disinfectant solution for the end user.

### COMPOSITION

Dimethyldioctylammonium chloride, Dipropylene glycol, 2-Aminoethanol, Tetrasodium ethylene diamine tetraacetate, Ethanol, Citric acid, Nonionic surfactant, Antifoam, Dye.

### INSTRUCTIONS FOR USE

Add 20ml ± 1 ml of concentrate to 980ml warm water (2% solution). Aspirate via the high and low volume suction tubes. Leave for a minimum of 15 minutes (to achieve full efficacy leave for 60 minutes) and then rinse thoroughly with clean water. Do not use on aluminium or copper parts. For professional use only.


### PACKAGING

#### Viroxid Aspirator

Code 070300..... (3) twin neck 1lt. bottle concentrate (2%)

All primary packaging are compliant to the technical specifications provided by the Pharmacopoeia. The material used does not contain latex and it is perfectly compatible with all the components of the formulation.

Medical Device Class IIb  
CE n° 1639



Not  
Foaming

Active  
in 15'

Active  
on Salmonella

### INFORMATION FOR THE EXCLUSIVE USE OF HEALTHCARE PROFESSIONALS

Manufactured by:  
Topdental (Products) Ltd.  
12 Ryefield Way, Silsden, West  
Yorkshire. BD20 0EF. England UK.  
Tel +44 (0)1535 652 750 www.topdental.co.uk

Distributed exclusively by:

**IDS S.p.A.**

Via S. Cristoforo 28/10


17100 Savona - Italy

www.ids dental.it - info@ids dental.it

### MICROBIOLOGICAL EFFICACY

MICRO-ORGANISM	CONC %	TIME
<b>Bactericidal</b>		
<i>Pseudomonas aeruginosa</i>	2%	15 mins
<i>Staphylococcus aureus</i>	2%	15 mins
<i>Enterococcus hirae</i>	2%	15 mins
<i>Legionella pneumophila</i>	2%	15 mins
<i>Salmonella typhimurium</i>	2%	15 mins
<b>Yeasticidal</b>		
<i>Candida albicans</i>	2%	15 mins
<b>Viruses</b>		
Hepatitis B	2%	15 mins
Hepatitis C	2%	15 mins
HIV	2%	15 mins
<b>Tuberculocidal</b>		
<i>Mycobacterium terrae</i>	2%	60 mins
<b>Fungicidal</b>		
<i>Aspergillus brasiliensis</i>	2%	30 mins

### PATHOGENIC AGENTS CLASSIFIED ACCORDING TO THEIR RESISTANCE TO DISINFECTANTS



### TEST STANDARDS

The following test standards have been achieved:  
EN 14885, EN 14476, EN 14348, EN 13727,  
EN 13697, EN 13624, ASTM E1052.

### PHYSICAL PROPERTIES

Yellow liquid with a characteristic odour  
Biodegradable  
pH: 11.0-12.0  
36 month shelf life from date of manufacture

### INFORMATION FOR THE EXCLUSIVE USE OF HEALTHCARE PROFESSIONALS

Manufactured by:  
Topdental (Products) Ltd.  
12 Ryefield Way, Silsden, West  
Yorkshire. BD20 0EF. England UK.  
Tel +44 (0)1535 652 750 [www.topdental.co.uk](http://www.topdental.co.uk)