ARCUS

EN

OPERATING INSTRUCTIONS MANUAL

LED Photoactivation Device

BLANCCONE[®] ARCUS[®] PRO

Ref.N 676200 (MOBILE) □

Ref.N 676201 (FIXED) □

BLANCCONE[®] ARCUS[®]

Ref.N 67600001N (MOBILE) □

Ref.N 67600002N (FIXED) □

Attention!
Before installing
and starting to work
with the device,
carefully read this manual!

I. DESCRIPTION AND FUNCTIONS.....	2	VIII. DAILY CARE AND MAINTENANCE.....	20
II. SYMBOLS	3	IX. PROBLEMS AND SOLUTIONS	20
III. SAFETY PRECAUTIONS.....	5	X. DISCLAIMER.....	21
IV. TECHNICAL DATA	9	XI. WARRANTY TERMS.....	22
V. COMPONENTS.....	10	XII. SERVICE DATA.....	23
VI. INSTALLATION	11	XIII. DECLARATION OF CONFORMITY	24
VII. DEVICE OPERATION	14	XIV. BLANCCONE ARCUS DATA	121

BLANCONE WHITENING SYSTEM

BlancOne ARCUS lamp and BlancOne photochemical whitening Gels (ULTRA, TOUCH, CLICK) are engineered to function as a cohesive system, essential for ensuring patient safety and treatment efficacy. BlancOne ARCUS lamp must be used exclusively with BlancOne photochemical whitening gels and not independently. Conversely, BlancOne photochemical gels should only be activated using BlancOne ARCUS lamp. For patient safety and to ensure the efficacy of the treatment, BlancOne photochemical gels (ULTRA, TOUCH, CLICK) must be activated using BlancOne ARCUS lamp. The unique properties of the Blanc One Arcus lamp, including specific wavelength and intensity, are essential for the proper activation of these gels. Patients must be protected using BlancOne CARE protections for eyes, face and lips.

Dental care professionals must undergo a training and certification to verify their understanding of and commitment to adhering to safety protocols when using the BlancOne In-Office Whitening System. Authorized dealer is responsible for providing or facilitating this training and certification, ensuring that all users of the

system are compliant and competent according local healthcare guidance.

Any serious adverse health issues or reactions associated with the BlancOne Whitening System must be reported to the authorized dealer who is responsible to record a comprehensive report about the incident and for reporting it to national healthcare authority and IDS Spa. Are exempted from this requirement common adverse effect of teeth whitening such as mild and transitory hypersensitivity and gum irritations.

BLANCONE ARCUS - LED PHOTOACTIVATION DEVICE

BLANCONE ARCUS - LED photoactivation device is a professional electronic device designed to emit specific wavelengths of light for use in professional dental settings. It provides a controlled light source within the 430-490 nm wavelength range.

For safety considerations, it is recommended to use eye, face, and lips protection during use. Dental professionals using the BLANCONE ARCUS should familiarize themselves with its proper operation and recommended safety guidelines. The authorized distributor is available to provide product information and user support as needed.

Any unexpected device-related issues should be reported to the authorized distributor, who will document the case and, if necessary, communicate with relevant national authorities. However, temporary effects commonly associated with whitening procedures, such as mild tooth sensitivity or mild and transitory gum irritation, do not require reporting.

I. DESCRIPTION AND FUNCTIONS

BLANCONE ARCUS photoactivation device is a professional electronic device designed to emit specific wavelengths of light that activate a photochemical reaction within BlancOne in-office gels. It is intended for use by licensed dental professionals in a clinical setting.

This device provides a controlled light source optimized for the photochemical properties of BlancOne-formulated gels, specifically within the 430-490 nm wavelength range. The whitening effect results solely from the interaction of the light with the whitening gel, according to the gel manufacturer's instructions.

BLANCONE ARCUS is not a medical device and does not itself perform any whitening treatment, alter tooth structure, or induce a therapeutic effect. The device solely emits light to activate a photochemical reaction and does not have an independent clinical or diagnostic purpose.

The selection and application of the whitening gel, as well as the use of any light exposure during the procedure, remains under the professional responsibility of the dental practitioner, who must follow the gel manufacturer's usage guidelines. The device manufacturer ensures that the BLANCONE ARCUS complies with applicable safety and quality standards.

BLANCONE ARCUS consists of light emitting head, stand, and power adapter. It is a CE marked electronic device, manufactured by BG LIGHT PRO Ltd in compliance with the following standards:

- EN ISO 9001:2015 Quality management systems
- EN 60601-1 Medical Electrical Equipment General Safety
- EN 60601-1-2 Electromagnetic Compatibility

II. SYMBOLS

European Conformity

Indicates conformity with local laws and regulations within the European Economic Area.

Caution

To indicate that caution is necessary when operating the device or control close to where the symbol is placed.

Refer to instruction manual/booklet

To signify that the instruction manual/booklet must be read.

Electronic instructions for use

To indicate on product or product packaging that relevant information for use of the product is available in electronic form rather than, or in addition to, printed paper form.

Dangerous voltage

To indicate hazards arising from dangerous voltage.

Hazardous light emission

To indicate hazards arising from light radiation.

Hazardous thermal effects

To indicate hazards arising from thermal effects.

Manufacturer

Indicates the device manufacturer.

Date of manufacture

Indicates the date when the device was manufactured.

Catalogue number

Indicates the manufacturer's catalogue number so that the device can be identified.

Serial number

Indicates the manufacturer's serial number so that a specific device can be identified.

Batch code

Indicates the manufacturer's batch code so that the batch or lot can be identified.


Applied part type B according to electric safety classification.

For indoor use in dry conditions.


Protection against solid objects larger than 12.5mm; no protection against water. (Ingress Protection rating)

For indoor use only

To identify electrical equipment designed primarily for indoor use.

Temperature limit

Indicates the temperature limits to which the device can be safely exposed.

Humidity limitation

Indicates the range of humidity to which the device can be safely exposed.

Keep away from sunlight

To indicate that transport package shall not be exposed to sunlight.

Keep away from rain

To indicate that the transport package shall be kept away from rain and in dry conditions.

Fragile; handle with care

To indicate that the contents of the transport package are fragile and the package shall be handled with care.

This way up

To indicate the correct upright position of the transport package.

Stacking limit by number

To indicate that the items shall not be vertically stacked beyond the specified number, either because of the nature of the transport packaging or because of the nature of the items themselves.

Waste Electrical and Electronic Equipment (WEEE)

According to Directive 2012/19/EU, this symbol indicates that the product should not be disposed of as urban waste at the end of its operating life. It must be sent to separate collection facilities for recovery and recycling.

III. SAFETY PRECAUTIONS

GENERAL WARNINGS:

To ensure safe use for both staff and patients, the following rules must be observed:

1. Do not allow unauthorized and untrained personnel to use the device to avoid risks.
2. Disconnect the device from the mains after completing the procedures.
3. Do not use or store the device in a dusty environment.
4. Do not expose the device to direct sunlight.
5. Do not spray disinfectant directly into the device - only wiping with a tissue soaked in disinfectant is acceptable.
6. Do not get the device, cables, or adapter wet or allow liquid drops to fall on the device, to avoid electric shock or damage to the device.
7. Store the device in a dry place, moisture can cause electric shock and damage.
8. The emitting head must not be covered during operation, the cooling openings must not be blocked, to avoid overheating the device and ignition.
9. Strong electromagnetic fields in the building can cause interference and malfunction of the device. If their source cannot be determined, change the location of the device and plug it into another outlet or another room, even in another building.
10. Before each patient, the output window of the emitting head and the spacer must be disinfected (using a disinfectant).
11. Before patients undergo whitening procedure, they must sign an Informed Consent, in which the possible side effects of the process are described in great detail. You can use a sample text from www.blancone.eu
12. It is recommended that the dentist attend an appropriate course "Teeth Whitening in a Dental Office"!
13. The device should not be used if some of its functions or parameters do not work or there are doubts about such, contact the manufacturer or service.
14. Opening and repairing the device may only be carried out by authorized service technicians from the manufacturer.
15. When replacing defective parts, only original BLANCONE ARCUS PRO parts must be used. The device or any of its parts must not be disassembled while it is connected to the network!

- **Fragile!** When transporting, using and storing the device, be careful! When transporting in an assembled state, be accompanied by two people, holding the moving parts. In case of unevenness on the floor (surface), the device must be carried by two people.

- According to Directive 2012/19/EEC, this symbol indicates that the product must not be disposed of as municipal waste at the end of its service life. The product must be taken to a specialized center for the separate collection of electrical and electronic equipment according to local legislation. Proper disposal of equipment that is no longer in use prevents negative consequences for the environment and human health!
- All packaging materials of the product must be kept out of reach of children to avoid risks of injury/suffocation.

RISKS

The device must be used in strict accordance with the Operating Instructions Manual.

1. Electrical safety

Before starting the device, make sure that the voltage and the type of plug correspond to the mains supply in the country. Use only the original adapter type FSP060-DAAN3 or SINPRO HPA65A-240.

Electrical safety is ensured by class I protection against electric shock according to EN 60601-1.

BLANCONE ARCUS must only be operated indoors, under the following conditions:

- temperature from + 10 ° to + 40 ° C;
- relative humidity 30 - 75%;
- lack of dust in the room;
- atmospheric pressure 700 - 1060 hPa;
- absence of chemically active and flammable substances;
- no part of the device should be wetted or immersed in water;
- the device or any of its parts must not be disassembled while it is connected to the mains!

To avoid the risk of electric shock, this device must only be connected to power supplies with protective earthing. Protect the device cables from damage to the insulation and tearing by sharp objects, strong pulling, rodents, chemicals, high temperature. In the presence of such damage to the electrical cables, the device must be repaired in the company service. The device must not be used with damaged cables.

During thunderstorms, it is mandatory to stop the procedures and disconnect the plug from the power supply.

Risk: Failure to comply with these instructions may result in electric shock to the device users.

2. Light radiation

BLANCONE ARCUS is a source of intense light in the blue range, to which the human eye has a high sensitivity. This results in taking serious measures for protection of patients, medical staff and random people nearby.

Mandatory use of protection goggles for the operator, and for the patient – goggles, BlancOne cheek retractor and lip protection (apply vaseline on lips for better comfort) in conjunction with BlancOne face protection paper.

Irradiation of the eyes and skin with intense light carries a risk of damage from light and heat. Freckles, pigmentation, redness and burns on the skin and gums are possible.

Never direct the light at the eyes! Irradiation must be limited to the work area.

It is mandatory to use special protective glasses from the set that meet the requirements of EN 166 and:

- tightly cover the eyes and temples, even if the person is wearing optical glasses.
- be made of volumetric colored impact-resistant plastic.
- do not transmit light with a wavelength of 430 - 490 nm.
- reduce the intensity of the blue spectrum by more than 100 times.
- have a stable mechanical structure, do not have scratches, cracks and damage on their surface.

The device can be used only after medical consultation on or by persons suffering from photo-biological reactions; persons taking photosensitive medications; persons operated on for cataracts, persons with retinal diseases, etc.

The risk of improper irradiation is severe eye irritation, temporary spots in the visual field, severe vision damage in case of direct irradiation, up to loss of vision.

3. Thermal effect

Thermal effect due to the absorption of blue light energy in tissues, in which the energy is converted into heat. The risk is only with prolonged overdose of irradiation.

Risk of pain, burning of soft tissues.

4. Fire safety

- Keep the device away from solvents, flammable liquids and heat sources.
- Do not expose the device to direct sunlight.
- Do not allow liquids and cleaning agents to penetrate the device, because this can lead to a short circuit and fire or cause functional failure.
- If the product emits an odor or smoke - disconnect it from the network, do not attempt to repair it, have it repaired at an authorized service center.

Risk of injury, burns, fire and damage.

5. Danger from mechanical moving parts

- Use in rooms with a horizontal floor surface.
- Mechanical parts should be assembled, positioned, and locked slowly and carefully.
- The device should not be used for transporting or moving people and objects.
- Always use two hands to position the stand in front of the patient's face.
- The radiation head should be fixed well - at an appropriate height and distance so that it is stable in front of the patient's face during the procedure.
- Take measures against damage to the human body from mechanical parts (moving and fixed), pinching, inertial reinforcement of the structure, manipulations with the device should be carried out carefully to avoid injury.
- Movement of the device parts without first unscrewing the fasteners may damage their mechanism. Before moving the device, the wheel stoppers must be released.

Risk of mechanical damage to the operator, the patient, and the device.

6. Medical and general prohibitions for using the device:

6.1 The device should not be used by:

Pregnant and nursing mothers, patients under 18 years of age, patients with severe periodontal pathology, with recessions, dental hyperesthesia; patients with allergies, wounds and infections, fresh scars on the face, skin infections, recently placed dental implants or surgical procedures in the oral cavity and face, fibrotic conditions, herpes, bleeding, bruises, burns, oncological diseases or indications for such on the face and oral cavity, atypical warts in the irradiation area, wounds that are difficult to heal; patients taking painkillers that dull the skin's sensitivity to heat; persons under the influence of alcohol or narcotic substances or those in an inadequate state.

6.2 The device can be used only after consulting a specialist doctor regarding the following cases:

Persons suffering from photobiological reactions; persons taking photosensitive medications; persons with allergies, persons who have recently undergone cosmetic surgery on the face or lips, including hyaluronic acid or Botox injections; persons with very sensitive skin, etc. If photosensitizing agents or medications are being taken, check the package insert of the medicine and do not undergo bleaching procedure if it is stated that it may cause photoallergic reactions, if after taking this medicine you are required to avoid sun exposure.

7. Procedures before using the device:

1. Conversation with the patient to clarify his status, whether he is suitable for the procedure and to familiarize himself with the possible risks and prohibitions, as well as to answer all questions regarding teeth whitening.
2. If the status result is positive, the content of the "Informed Consent" is continued and the patient must sign it.
3. Hygiene procedures - cleaning of tartar and, in a short time, arranging an appointment for the whitening procedure.
4. It begins with isolating the soft tissues with a retractor, a protective face napkin, protective cream/vaseline on the lips.
5. Mandatory use of the protective glasses provided by the manufacturer for the operator and the patient. Do not use other types of glasses, because they may not provide the necessary protection for the eyes.
6. Application of a gingival barrier, polymerization and application of the whitening system/gel according to its protocol.
7. Do not irradiate the patient by contact, but from a distance determined by the spacer.
8. Check the patient's sensitivity to the light intensity and select an appropriate operating mode.
9. It is important to check how the patient reacts to this intensity and adjust it so that there is no risk of an adverse reaction.
10. There should be continuous control over the procedure to ensure that the distance has not changed. The patient should not be left unattended during the procedure. If the patient is uncomfortable, the procedure should be stopped for a while or terminated.
11. It is desirable to have a photo before and after the procedure of the teeth and soft tissues and to monitor whether there is redness or other reaction. The photo should be stored for a long enough time to monitor the whitening effect and possibly the reaction of the soft tissues.

8. Procedures after using the device:

1. Examination of the patient after the procedure for redness and changes in the mucous membrane. If there is any, give a prescription for appropriate treatment and maintain contact with the patient in the following days until the possible problem disappears. It is desirable for the dentist to maintain contact with a dermatologist and to be able to offer a consultation with one if problems arise with the patient.
2. The patient should be warned to observe the necessary hygiene and not to undertake self-treatment, which may worsen the reaction and to maintain contact with the dentist if necessary. Do not eat aggressive foods with high acidity and pungency.
3. The clinic / dentist / takes full responsibility for quickly and effectively solving any possible problems in order to obtain the desired aesthetic result without unwanted complications and consequences.

IV. TECHNICAL DATA


1. Operating voltage to the power adapter - 100-240V / 50-60Hz,
to the light emitting head - 24VDC.
2. Charging adapter type FSP060-DAAN3 or SINPRO HPA65A-240
3. Dimensions:
 - Light emitting head - length - 220 mm, width - 155 mm, height - 95 mm
 - Stand (mobile) - length (assembled) - 1270 mm
 - Stand (fixed) - length (assembled) - 52 mm / 125 cm (arm extended)
4. Weight:
 - Light emitting head - 600 g
 - Stand (mobile) – 8000 g
 - Stand (fixed) – 1700 g (without bracket)
5. Irradiance - light intensity:
 - BLANCONE ARCUS PRO - 100% / 75% / 50% (Custom mode)
 - BLANCONE ARCUS+ - HI (100%), LO (50%)
6. Operation time - 1 - 20 min / $\pm 5\%$ /
7. BLANCONE ARCUS PRO: Pre-programmed modes for CLICK, TOUCH & ULTRA whitening systems, for Fast / Standard / Extended / Sensitive whitening procedure
8. Emitted light - blue 430 - 490 nm.
9. Light source - 10 LED modules x 5W = 50W
10. Fan air cooling
11. Ability to start and stop the light emitting head at any time.
12. Degree of protection against electric current - applied part type B

The manufacturer of this device will provide upon request additional technical documentation / information necessary for the technical staff of the user to repair those parts of the device that the manufacturer has specified as subjects for repair.

V. COMPONENTS

MOBILE

1. Power cord
2. Power Adapter
- 100-240VAC / 24VDC
3. Straight & curved arm
4. Back fixator
5. Front fixator
6. Base with 5 wheels
7. Vertical arm
8. Light emitting head
9. Distance stick (spacer)
10. Control panel
- ARCUS PRO / ARCUS+**
- 10A. START/STOP/OK button
- 10B. PAUSE / RETURN button
- 10C. SCROLL DOWN button
- 10D. Display
11. Emitting window of the light emitting head
12. Cables and connectors
13. Safety glasses - 2 pcs


The device should only be moved when all 5 wheels are attached to the Base.

Each wheel has a locking mechanism.

Positioning of the horizontal arm is done by the back fixator.

FIXED

1. Power cord 24VDC to head
2. Straight arm
3. Curved arm
4. Front fixator
5. Back fixator
6. Fixing bracket
7. Emitting head
8. Distance stick (spacer)
9. Control panel
- ARCUS PRO / ARCUS+**
- 9A. START/STOP/OK button
- 9B. PAUSE / RETURN button
- 9C. SCROLL DOWN button
- 9D. Display
10. Emitting window of the light emitting head
11. Cables and connectors
12. Power adapter with mains cable
13. Vertical pole of dental unit
14. Safety glasses - 2 pcs

Fine positioning is done by turning the light emitting head so that the light is perpendicular to the patient's teeth and centered on them.

If the movement of the light emitting head is difficult or it does not hold the selected position it is necessary to adjust the Front fixator.

VI. INSTALLATION

MOBILE

BLANCONE ARCUS PRO / ARCUS+ is intended for use as a stand-alone device. Use the device indoors and only on a horizontal surface.

1. Unpack the transport boxes, remove the components and all their protective covers.

2. Assemble the base with five wheels.

3. Mount the vertical arm to the center of the base.


4. Insert the cable coming from the curved arm into the top opening of the vertical arm until the coupling appears at the bottom of the wheelbase.

5. Mount the curved and straight arms to the vertical arm and fix it with the side stop screw using a 3 mm Allen key (included in the package).

6. Tighten the side stop screw on the wheeled base using 3mm Allen key to fix the vertical arm into the base.
7. Disassemble the cover of the power adaptor box located on the vertical arm by gently lifting it upwards and pull it.
8. Connect the coupling coming from the bottom of the base with the power adapter and connect the power cord through the opening on top of the box.
9. Insert the power adapter in the designated box and fix it with the cover

FIXED

10. Unpack the transport box, remove the components and initially assemble the fixing bracket to the dental unit. Make sure the diameter of the bracket matches the diameter of the dental unit vertical arm.

Soft PVC tape strips are provided in the case that there is a need to compensate for small difference in the two diameters. Tighten the two screws with a 4 mm Allen key (included in the package).

11. Tighten the stop/safety screw with a 3 mm Allen key (included in the package).


12. Insert the cable through the small opening of the fixing bracket and place both arms in it.

MOBILE / FIXED

13. Tighten both front and back fixators.

14. With the front fixator firmly tightened, mount the light emitting head.

15. Connect the coupling coming from the light emitting head to the one from the curved arm.


16. Connect the cable on the bottom of the curved arm with the power adapter using the connectors.

17. Plug the power cord into the power network.


Always ensure correct connection!

The device will not work if the cable connectors are incorrectly coupled.



VII. DEVICE OPERATION

1. The device is placed in a suitable position for operation:

2. The emitting head is aligned with the dental arch:

3. Use the spacer located on the underside of the emitting head to measure the correct distance to the patient's teeth.

EMITTING HEAD – BLANCCONE ARCUS PRO

BUTTON	FUNCTION
	ARCUS PRO START / STOP – long press – 2-3 sec. WAKE UP – long press – 2-3 sec. OK – confirm the operation (short press)
	SCROLL DOWN – navigate the menu from top to bottom (short press) QUICK SCROLL – when selecting time in Custom mode (long press) ADJUST INTENSITY – while working in Custom mode (during active work cycle)
	PAUSE – temporarily stop operation; restart after that RETURN / BACK – return to the previous menu

Starting the device

When turning on the display shows the logo and a melody is heard.

A reminder screen to sign the


Patient Informed Consent

(link with QR code is available in menu INFO > Informed consent)


A reminder screen to protect the patient:

safety glasses, retractor, protective face napkin

Selection of bleaching gel **CLICK / TOUCH / ULTRA** , gingival barrier or custom mode


The options are scrolled with the “arrow” button,
confirmation of the gel selection is made with
the “Start/Stop” button.


Selection of whitening procedure: **Fast / Standard / Extended / Sensitive**

Example:

ULTRA⁺ gel / Standard mode is selected

Graphic of selected program


Active cycle from selected program

Prompt screen to **reapply gel** between cycles.

During this screen, the device is in a
“technological pause”.

Switching to the next mode is done by pressing
the “Pause” button.

After the 3 cycles have elapsed,
the program is finished.


At the end - screen that **all cycles are finished**

In the same way, you can choose **CLICK**, **TOUCH** gel and one of 4 whitening procedure programs (**FAST**, **STANDARD**, **EXTENDED**, **SENSITIVE**).

Polymerization of **gingival barrier**


(fixed 1 cycle – 1 min,
at 100% light intensity)

Custom settings

Choice of light intensity (50 / 75 / 100%)


Choice of timer (1 – 20 min)

Example: active cycle: 8 min, 75%


Info menu

If it is necessary to change the language, this can be done from the INFO > LANGUAGE menu.

If **overheating** occurs during device operation (thermal protection is activated), the device stops light emitting and a screen appears:

The screensaver is activated after 15 minutes of inactivity. Pressing any button returns to the last selected operating screen. Turning off the device – by long pressing the START/STOP button.

After that disconnect the device from the mains.

A recommended check of the light intensity is required after a number of operating hours of the LED modules specified by the manufacturer. The display screen will show warning message: "Light intensity check recommended - contact your local dealer"


EMITTING HEAD – BLANCONE ARCUS*

ARCUS*

1. The device is placed in a suitable position around the patient. Emitting head must be in the same plane as the dental arch. The appropriate distance to the patient's teeth is determined using the spacer.
2. Plug the power adapter into the power network. Only one dot on the display lights up. Press START / STOP button and then the "Hi" (High) or "Lo" (Low) intensity indication appears.

3. The display shows "HI" - high intensity level. If you want to switch to the "Lo" (Low) intensity mode, press the TIMER / MODE (+ / -) button once. Every time before setting the required time, you will be able to select the Hi or Lo intensity level. Press START / STOP button to confirm your choice.
4. Select timer by pressing TIMER + (increasing time) and TIMER- (reducing time) - between 10 sec. – 20 min. The display will indicate timeframes between 10 sec. and 9 min. and 50 sec. with a dot between the two digits. The device will save the time frame used in the last procedure.
Examples:
Display indicates: 13. Timer is set to 13 min.
Display indicates: 0.1. Timer is set to 10 sec.
Display indicates: 9.5. Timer is set to 9 min. and 50 sec.
5. Press the START / STOP button one final time and the device will start to emit light. A flashing dot on the display will indicate that the device is working properly.
6. If any of the TIMER / MODE (+ / -) buttons is pressed during operation, the unit switches to Pause mode. A blinking display with the remaining time will indicate that the device is paused. Resume operation by pressing any TIMER / MODE (+ / -) button again.
7. To Stop the device at any time press the START / STOP button. If the device is not manually stopped, it will stop after the selected time has elapsed. The cooling fan will continue working for some time.
8. At the end of the day turn off the device by holding down the Start / Stop button for 3-4 seconds. After that unplug the adapter from the mains.

9. In case of overheating the thermal protection is triggered and the device will stop emitting light. The display will indicate "Oh" and the fan will turn on for 1 minute. After cooling, the device can continue its normal operation, staying in Pause mode.

VIII. DAILY CARE AND MAINTENANCE

1. Disinfection of the PVC window of the radiation head

For each patient, the radiation window is cleaned with cotton wool slightly soaked in an alcohol solution.

2. Cleaning the device

To disinfect the device and its parts, spray disinfectant on a piece of soft cloth or cotton wool and clean. Do not spray directly on the device!

Do not use abrasives or solvents, they will damage the device parts!

3. Cleaning the spacer

The spacer can be easily rotated, disassembled (by pulling) and cleaned with a disinfectant/ alcohol solution.

Do not disinfect the spacer by spraying directly on it (without disassembling it).

IX. PROBLEMS AND SOLUTIONS

The device does not work

Check whether the device is plugged into the outlet and whether the power cord and connectors are connected correctly. Use only the original adapter type FSP060-DAAN3 or SINPRO HPA65A-240. If the adapter has been dropped, mechanically damaged or wet, it may not work. Replace with a new one!

The device blinks and does not work

With a damaged adapter, the device may work properly until the blue light starts. Then it blinks for a moment and stops, and it is not possible to perform whitening. Reason could be dropped, mechanical-

ly damaged or wetted adapter. Replace the adapter with a new one.

Damaged power cord

Do not use the device. Turn off the power and replace the cord with a new one.

The emitting head cannot be fixed

Tighten the fixators once the desired head position is determined.

When starting the device, the display shows E1 (ERROR 1, or other) and goes out

The device should be sent to the manufacturer's service center.

X. DISCLAIMER

The manufacturer BG LIGHT PRO Ltd is not responsible and does not owe compensation for damage or lost profits in the following cases:

1. If the device is not installed correctly or its components are not used for the intended purpose and method of operation by the manufacturer, which causes harm to the patient, the dentist or another person.
2. The expected effect of the procedure is not achieved.
3. In case of injury to the patient, for example – burning and swelling of soft tissues, swelling of the lips, appearance of freckles, pigmentation, etc., which may occur if the requirements are not followed (points 4-15 of this section).
4. The patient is not familiar with and has not signed an "Informed Consent" before the procedure, which includes "Risks" p.6.1, p.6.2 of Chapter III.2 of the Operating instructions manual of the device. It can be downloaded from our website - www.blancone.eu
5. If the patient is not prepared for the procedure. At least 24 hours have not passed since the tartar removal. Soft tissues have not been isolated with a gingival barrier, retractor and face napkin!
6. If the whitening system has come into contact with unprotected areas of the gums, mucous membrane, tongue, lips, face. Lack of photographic material of the patient before starting the whitening, which would certify that the protocol for preparing the patient before whitening, described in Chapter III, point 7 of the Operating instructions manual of the device, has been followed.
7. If the dentist has not followed "Safety Measures and Risks" p.6, p.7, p.8 of Chapter III of the Operating instructions manual of the device.

8. If the dentist left the patient unattended during the procedure and despite symptoms of discomfort and pain, did not terminate the procedure.
9. In case of radiation overdose (see Chapter VII - types of programs).
10. Use of the whitening system / gel that does not require light activation.
11. Use of the whitening system that is activated with a wavelength other than 430-490 nanometers.
12. Use of the whitening system with an expired expiration date.
13. When the procedure is performed by untrained personnel.
14. When treating facial and dental tissues with substances that would cause photosensitivity or allergic reactions during light photoactivation.
15. If the protection napkins, retractors, aids and instruments used are treated with disinfectants or other chemical preparations that would cause allergic reactions during photoactivation.
16. The manufacturer reserves the right to change the product in order to improve quality and functionality.

Note: Please check the website www.blancone.eu for the latest version of the “Instructions for Use” and “Patient Informed Consent”. We are committed to maintaining this information according to the latest developments in this field.

XI. WARRANTY TERMS

1. The warranty period of the BLANCONCONE PRO whitening activation device is 24 (twenty-four) months from the date of purchase. It is verified with a complete warranty card or invoice. If these documents are not provided, the date of manufacture is accepted as start of the warranty period.
2. During the warranty period, the replacement of defective elements is carried out free of charge by the manufacturer.
3. If during the warranty period the device becomes defective due to improper operation (mechanical, chemical, thermal, electrical damage, improper storage) and other reasons due to the user's fault, the repair is not covered by warranty.
4. If liquids such as water and solvents, aggressive and flammable substances and their vapors, penetration of foreign bodies or others into the device have entered the device, it must be disconnected from the power supply and not used. The company service will determine the damage, as the repair is not covered by warranty.

5. No claims for damage and complaints are accepted as a result of electric shocks, thunderstorms or failure to comply with electrical safety measures.
6. When repairs are carried out by unauthorized people outside the company service or if non-original elements and assemblies are used, the user loses the right to a free warranty service.
7. Sending for repair to the company service must be in the original packaging of the product and is at the expense of the sender.
8. Warranty service is carried out at the company service at the address:

IDS SPA

Via San Cristoforo, 28/10 - 17100 Savona

Tel. 019-862080 - www.ids dental.it

9. The device's service life is set at 5 years from the date of purchase. After the end of its life cycle, the user must discontinue use of the device.

At the user's request, the life cycle can be extended. For this purpose, the device must be returned to the manufacturer BG LIGHT PRO Ltd, where its characteristics must be checked and all main elements (emitters, adapter, electronic boards) tested and replaced.

10. In the event of disputes arising in connection with the application and interpretation of these Operating Instructions, they will be resolved by the Plovdiv District Court, in accordance with the current Bulgarian legislation.

XII. SERVICE DATA

Repair service only valid for the EU territory, for countries outside the EU please contact your dealer. All the equipment distributed by us that will need technical interventions during or outside the warranty period, must be sent directly to the company Vi.Vi.Med srl who will manage the repair and invoice.

VI.VI.MED Srl **Technical Service**

Via Tevere 11 - 51018 Pieve a Nievole (PT) – Italy

Tel. +39 0572 520910 - fax. +39 0572 520807 - Mobile: +39 366 6362021

e-mail: info@vivimedsrl.com – web www.vivimedsrl.com

XIII. DECLARATION OF CONFORMITY

 TECHNICAL FILE BLANCONE ARCUS LED photoactivation unit	EU Declaration of conformity	TD 5.3
		<i>Revision 03</i>

Manufacturer: **BG LIGHT PRO Ltd**

Address: 155, Vasil Aprilov blvd., 4027 Plovdiv, Bulgaria
 Tel.: +359 32 644089, +359 888 809256, email: office@bglight.com
 BULSTAT UIC 208000520, VAT N: BG208000520

Product:	Product code:	Name:
LED photoactivation unit	676200 600-003(A)IT	BLANCONE ARCUS PRO (MOBILE)
	676201 600-003(BT)IT	BLANCONE ARCUS PRO (FIXED)
	67600001N 600-001+	BLANCONE ARCUS+ (MOBILE)
	67600002N 600-001bt+	BLANCONE ARCUS+ (FIXED)

Intended purpose: BLANCONE ARCUS is designed for photoactivation (irradiation of blue light 430-490 nm) of whitening gel based on hydrogen or carbamide peroxide, intended for light activation, used in dental practice.

The manufacturer declares under his own responsibility that the described product complies with the essential requirements of the following normative acts and regulatory technical documents when used for its intended purpose and in accordance with the safety instructions:

Document	Title	Edition / date of issue
<i>Directive 2014/35/EC</i>	<i>Electrical equipment designed for use within certain voltage limits</i>	26.02.2014
<i>Directive 2014/30/EC</i>	<i>Electromagnetic compatibility (EMC)</i>	26.02.2014
<i>Directive 2011/65/EU</i>	<i>Restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)</i>	08.06.2011

To achieve compliance, the requirements of the following harmonized standards are met:

Harmonized standard	Title	Edition / date of issue
EN 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance. IEC 60601-1:2005/ IEC 60601-1:2005/A1:2012	2006+ A1:2013;A12:2014;A2:2022
EN 60601-1-2	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance. IEC 60601-1-2:2007	2015+ A1:2021
EN 60601-1-6	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability. IEC 60601-1-6:2010	2010+ A1:2015; A2:2021
EN IEC 63000	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances (IEC 63000:2016)	2019
EN ISO 9001	Quality management systems - Requirements	2015
EN ISO 14971	Medical devices – Application of risk management to medical devices. (ISO 14971:2019)	2019+ A1:2022


The products described above comply with the essential requirements specified in Directives 2014/35/EU, 2014/30/EU and 2011/65/EU.

The declaration of conformity is issued in fulfillment of the requirements for the declaration of conformity, according to the relevant annexes of the implemented directives, based on test results - protocols No. 24.0027/02.057 - 27.01.2025 and 2EMC-C-23-041 / 17.01.2025, the assessed conformity and the maintained Internal production control system - Quality management system - Module A - ISO 9001:2015 certificate – No AC090 100/2409/6029/2025.

BG LIGHT PRO Ltd maintains a Technical File according to the requirements of the described directives.

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BG LIGHT PRO Ltd posee un archivo técnico de acuerdo con los requisitos de las directivas descritas.