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Harvard Temp C&B is a composite material paste, ready-to-use in double-chamber cartridge and based on methacrylates. It allows the confection of temporary crowns and bridges, inlays, onlays and veneers, also directly on the patient.

Thanks to its high final hardness and to its elasticity, the material is particularly suitable for confecting multiple bridges. Cured material demonstrates a very good stability, especially at the thin edges of the crowns.

Harvard TEMP C&B is free of methyl methacrylate. Its curing temperature stands under 40 °C / 104 °F.

After seating of the provisory, prepared teeth are protected against extern attacks and occlusion is maintained. The cured provisory shows very good values of flexural strength and abrasion resistance, a low polymerization shrinkage and a high fitting accuracy. It is also characterized by a good cosmetic result, an optimal polishing, high color stability and natural fluorescence.

Failures by confecting the provisory with Harvard TEMP C&B can be easily corrected, inasmuch as already hardened material aggregates freshly extruded material and light-curing composites.

Composition

Unsaturated esters, multifunctional methacrylates, multifunctional acrylates and malonyl urea.

Contra indication

Do not use with patients who show an allergy to any of the ingredients of Harvard TEMP C&B.

Processing

1. Preparation of the situation impression
Confect a situation impression using either impression materials from addition cross-linked silicones, or alginates.
For a better stability of the provisory, carve out or cut out the interdental areas.
2. Preparation of the cartridge
Remove the cap of the Harvard TEMP C&B cartridge and throw it away. Replace it by the 10:1 mixing cannula. Turn the cannula 90° until it locks in position.
Load the Dispenser (application gun) with the cartridge.
Remark: discard the first 2-3 mm of the extruded material. This must be done before each new mix. The used cannula plays the role of the cap until the next mix.
3. Application
 - 3.1 Load the situation impression with Harvard TEMP C&B to the appropriate positions. Extrusion occurs by slight pressure of the gun trigger. Plunge the mixing cannula into the material, in order to avoid the formation of air bubbles.
 - 3.2 Place the impression in the mouth.
Steps 3.1 and 3.2 must be done within 50 seconds (working time).
 - 3.3 After 1-2 minutes of setting time in the mouth, remove Harvard TEMP C&B together with the situation impression. In case of setting using a model, remove after 3-4 minutes.
The provisory should be at this step in a plastic state, i.e. cured but still elastic.
4. Post curing and finishing
After removing of material excesses, the provisory is post cured in warm water at 45–55 °C / 113–131 °F for a few minutes and can be then worked out.
In case of a post curing at room temperature, the provisory can be worked out and polished with rotary instruments only after a curing time of 6 minutes.
Do not breathe polishing dusts! Use a suitable mouth and nose protective device and/or work with an aspiration system!
The inhibition layer, which develops onto the surface of the provisory in presence of the atmospheric oxygen, should be removed with a suitable solvent, e.g. ethyl alcohol, before working out.
5. Luting
Luting of the provisory should be made using an eugenol-free temporary cement, e.g. Harvard TEMP, in order to avoid any influences on the curing of resin-based fillings and/or luting materials, which will be applied later.
6. Repairing
If, against all odds, a fresh-made provisory comes to a fracture, the broken parts can be coated and agglutinated with freshly extruded Harvard TEMP C&B.
In the case of a longer existing provisory, clean and roughen the areas of fracture, in order to provide them with mechanical retentions. Then join the fractured areas together after coating with freshly extruded Harvard TEMP C&B. Curing time can be accelerated in warm water (cf. § 4).

Process overview:

Up to 50 s	impression loading at 23 °C / 73 °F and application in mouth
50 s – 3 min	setting in mouth at 37 °C / 99 °F
3 – 4,5 min	polymerisation at ca. 50 °C / 122 °F
From 4,5 min	trimming / polishing at 23 °C / 73 °F

Warning advice

- Avoid contact with skin, mucous membrane and eyes!
- In case of contact with skin, immediately wash with water and soap!
- In case of contact with eyes, immediately rinse with a lot of water and seek medical advice!
- A sensitization to Harvard TEMP C&B may occur with some patients.
Do not use Harvard TEMP C&B anymore if allergic reactions are observed.
Do not use Harvard TEMP C&B with patients who show allergic reactions against acrylates.

For professional use only!

- Ordering-no.:
- 7081652 50 ml cartridge with 10 mixing cannulas 10:1, shade A2
 - 7081653 50 ml cartridge with 10 mixing cannulas 10:1, shade A3
 - 7081654 50 ml cartridge with 10 mixing cannulas 10:1, shade A3,5
 - 7094000 50 mixing cannulas 10:1
 - 7095000 Dispenser (application gun)

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