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## EN OPTICORE N3 Dual Curing Flowable Hybrid Composite for Core Build-ups and Cementation of Posts

Product description
OPTICORE N3 is a dual cure high radiopaque flowable hybride composite material with nano particles (nano-zirconium dioxide, nano-calcium-fluoride) for core-build-up and cementation of posts. Due to its excellent mechanical properties final crown preparation can be carried out more precise. OPTICORE N3 is based on poly- and difunctional methacrylates and inorganic filler particles of 0.02-10 µm. The total filler content is 64 % by weight and 48 % by volume. Delivered in 1:1 MINIMIX-syringes and AUTOMIX-cartridges, OPTICORE N3 can be easily dispensed and applied directly. OPTICORE N3 exhibits a short setting time without high heat generation. The dual cure properties enables the dentist, also to carry out cementations and core build-ups in cases where a light cure cannot be guaranteed to be sufficient.

OPTICORE N3 meets the requirements of DIN EN ISO 4049, type 2, class 3.

Indications/Intended use

- Core build-ups
- Cementing of posts
- Cementing of crowns and bridges, inlays and onlays

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

The placement of OPTICORE N3 is contraindicated if a dry working area or the recommended application technique are not possible.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group

Persons who are treated during a dental procedure.

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OPTICORE N3 may be contoured by using a composite instrument. A flat-ended interproximal carver is recommended.

Place OPTICORE N3 directly into the preparation and allow the system to self cure for 3:30 minutes. After that the material should be light-cured (40 seconds) with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm<sup>2</sup>. With this technique an optimum of physical properties will be obtained. An explorer can be used to test that the OPTICORE N3 is completely set. Remove the matrix not earlier than the material has set. Final core preparation preparation on the OPTICORE N3 can be carried out by using crown preparation burs. (see diagram 9)

**Additional Notes/Warnings**

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of the composite.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Color stability meets the requirements of DIN EN ISO 4049.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

Composition

Dimethacrylates, dental glass, silicone dioxide, photo initiators, catalysts

**Storage**

Do not store above 20 °C (68 °F). Protect from direct sunlight.Store unopened material in the refrigerator.

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Opened cartridges have to be used up within 3 months. Do not use after expiry date.

**Disposal**

Disposal of the product according to local authority regulations.

**Reporting obligation**

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

**Note**

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – https://ec.europa.eu/tools/eudamed).

**Warranty**

SanaPro Dental GmbH warrants this product will be free from defects in material and manufacture. SanaPro Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and SanaPro Dental GmbH's sole obligation shall be repair or replacement of the product.

**Limitation of Liability**

Except where prohibited by law, SanaPro Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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