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## GB) RECOMMENDATION FOR USE:

Core build-ups and cementation of posts with Opticore System

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#### 1) PREPARATION

adjust the post cavity depth (approximately 5-7 mm) with corresponding precision drill.

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#### 2) ETCHING

Etch dentine surfaces with phosphoric acid for 15 seconds. Wash and dry

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#### 3) DUAL BONDING

- place 2 layers of Opticore Dual Bond mixed at 1:1 ratio with catalyst. Dry with oil free air.

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#### 4) CEMENTATION

Place OPTICORE directly into the cavity using the Oral Tip with End Tip.

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#### 5) ENDODONTIC POST PLACEMENT

Insert carbon or glass post applying enough strength. Light cure for 20 seconds in order to hold the post in position.

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#### 6) CORE BUILD-UP

Build up core with Opticore. Light cure for 40 seconds or leave self-curing for 3:30 minutes.

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**COMPARISON TABLE:** This table shows that OPTICORE is much more flexible compared to competitor products, but similar hardness and compression strength, thus allowing a better overall strength to bite forces into the mouth. Also, thanks to its lower radio opacity, ensures a better esthetic performance under metal free crowns.

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# OPTICORE Classic Dual Curing Flowable Hybride Composite for Core Build-ups and Cementation of Posts

Rev3 - 27 Jan 2009

OPTICORE is a dual curing radiopaque flowable hybrid composite based material for core-build-up and cementation of posts. Delivered in auto-mixing 1:1 cartridges it can be easily dispensed with an application gun and applied directly. OPTICORE exhibits a short setting time without high heat generation. The dual cure properties enables the dentist, also to carry out cementations and core build ups in cases where a light cure cannot be guaranteed to be sufficient.

OPTICORE meets the requirements of the international standard ISO 4049.

## Preparing of the cartridge

Take off the turn cap of the double cartridge and throw it away (do not use it again!). It is substituted by the special 1:1 mixing cannula by attaching it with a sideways turning movement, locking the cannula in place.

The application gun is easy to load with the prepared cartridge and is ready for application.

*Note: The initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix.*

Store used cartridge with fixed used mixing cannula in the dark.

**The working time (23°C (74°F)) of OPTICORE in the self cure mode is 1:30 minutes from start of mixing.**

## 1. Core-Build-Up

### 1.1 Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

### 1.2. Cavity Preparation

Remove all existing old restorations and decay from the tooth. If necessary place any pins or posts. Refer to directions of the selected post manufacturer.

### 1.3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of hard setting calcium-hydroxide material.

### 1.4. Etching

Recommended is the total etch technique. Apply OPTIMA BLUETCH onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for 20 seconds. On primary teeth and teeth high in fluoride, a 60 second etch is recommended. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of OPTICORE DUAL BOND ADHESIVE. The etched enamel bevel should have a chalky white appearance.

### Etching Precaution:

*It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.*

## 1.5. Application of OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR

Apply a mix of OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR in a mixing ratio of 1:1 generously with a brush onto the etched enamel and dentine surfaces for 30 sec with agitation. Dry cautiously with oilfree air for 15 sec.

If root canals are to be primed - these are mostly not accessible by a brush - use a syringe for the application of OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR and flush it into the prepared root canals. The flushing must go on over a time of **60 seconds** discontinuous in **4 cycles**. All the dentin surfaces must kept wet with the primer over the specified time by generously brushing the accessible surfaces with a brush.

An excess of the primer in the canals has to remove with soft paper tips. Dry the root canal cautiously with oilfree air for 15 seconds and cure the coating with a dental halogen light unit for 20 seconds.

An application of a second layer of the dual cure priming and bonding mixture is strongly recommended before application of OPTICORE. For **30 sec** coat again all areas to be primed by the same way as described above. The dual cure system will autocure after application under the restorative placement within ca. **3 minutes** in cases where light cure cannot be guaranteed to be sufficient.

It is essential that the primed dentin and enamel surfaces are dry and contaminant free for the application of OPTICORE.

### 1.6. Application

Place the mixing cannula directly into the preparation and press out the paste.

OPTICORE is automatically mixed when dispensed with slight and even pressure. Filling should occur from bottom upward to prevent voids. To facilitate placement of OPTICORE place a matrix band around the prepared tooth.

OPTICORE may be contoured by using a composite instrument. A flat ended interproximal carver is recommended.

Place OPTICORE directly into the preparation and allow the system to self cure for 3 minutes. After that the material should be light-cured (40 seconds). With this technique an optimum of physical properties will be obtained.

An explorer can be used to test that the OPTICORE is completely set. Remove the matrix not earlier than the material has set.

Final crown preparation on the OPTICORE can be carried out by using crown preparation burs.

## 2. Post Cementation

### 2.1. Isolation

Use of a rubber dam to isolate the tooth is strongly recommended.

### 2.2. Root Canal Preparation

Refer to directions of the selected post manufacturer.

Remove all existing old restorations and decay from the tooth. Prepare and clean the root canal with hydrogen peroxide (3%) or sodium hypochlorite solution and remove excess solution from the canal with a soft paper tip.

Etch the root canal with OPTIMA BLUETCH for 15 seconds and rinse with water using an endodontic irrigation syringe. Blot the canal dry with soft paper tips, leaving the dentin visibly moist. A slightly wet dentine surface is important for the function of OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR

### Etching Precaution:

*It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.*

Flush the prepared root canal with the OPTICORE DUAL BOND ADHESIVE / OPTICORE DUAL BOND ACTIVATOR mixture by using a microapplication brush or a syringe. The flushing must go on over a time of **60 seconds** discontinuous in **4 cycles**. All the dentin surfaces must kept wet with the primer over the specified time by generously brushing the accessible surfaces with a brush.

An excess of the primer in the canals has to be removed with soft paper tips. Dry the root canal cautiously with oilfree air for 15 seconds and cure the coating with a dental halogen light unit for **20 seconds**.

### 2.3. Post Cementation

The OPTICORE DUAL BOND ADHESIVE/OPTICORE DUAL BOND ACTIVATOR mixture is applied generously with a brush onto the post. Dry cautiously with oil free air for 15 seconds. The prepared root canal is filled with OPTICORE by using a lentulo spirale or a paste filler. Seat the post into the canal gently and maintain firm pressure until the post is seated. OPTICORE self-cures within **3 minutes**.

As soon as the OPTICORE has set proceed with the core-build-up procedure.

## Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.
- Contact of resin pastes with skin and gingival tissue should be avoided, especially by anyone having known resin allergies.

## Storage

Do not store above 20 °C (68 °F). Store unopened material in the refrigerator. Opened cartridges have been used up within 3 months. Do not use after expiry date.