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Pre operator gel for dental hygiene and periodontal pocket cleansing

GB

INSTRUCTIONS FOR USE

Read all this leaflet carefully before you start using this product. Keep this leaflet, you may need to read it again.

1. WHAT IS ORALSAN OXYGEL AND WHEN USE IT

ORALSAN OxyGel is a pre operator gel based on 11% Carbamide Peroxide and 1% Disodium EDTA to be used sub and supra gingivally to dissolve plaque and soften tartar deposits.

Before starting the scaling procedures, the gel is easily flowed inside the sulcus by just placing the tip on the gum border and injecting the gel. In few seconds the gel will start releasing Hydrogen Peroxide which drastically controls the bio-load of the sulcus, while the combined action of Urea and Disodium EDTA will produce soluble complexes with calcium and heavy metal facilitating the dissolution of microbial plaque and tartar. This pretreatment will allow speed-up a lot the hygiene procedure working both with hand instrumentation and with no or minimal water spray during ultrasonic scaling, thus dramatically reducing the production of aerosols.

Hydrogen peroxide has a mild hemostatic action. ORALSAN OxyGel formula also includes Sodium Hyaluronate to improve the healing of gums after prophylaxis and Potassium Nitrate to reduce the sensitivity from the opening of dental tubules during the scaling. Dental hygiene session will be much faster and comfortable for the patient.

ORALSAN OxyGel can be used to brush teeth surface with a prophy brush or cup to remove effectively biofilm and prepare the tooth surface. This will speed up stain removal with air prophy and improve active oxygen penetration during in office teeth whitening procedures.

2. HOW TO USE

Connect the application tip with syringe and apply a small quantity of the gel into the periodontal pocket and on tooth surface. leave the gel to act for 1-2mins and after that start the cleansing procedure working manually or with minimal or no water spray to avoid the production of aerosols.

ORALSAN OxyGel can be used to brush teeth surface with a prophy brush or cup to effectively remove biofilm and prepare the tooth surface. This will speed up stain removal with air prophy and improve active oxygen penetration during chairside teeth whitening procedures.

3. POSSIBLE SIDE EFFECTS

If you experience any of the following serious effects STOP using this product immediately: Urticaria, skin rash and itching, allergic reactions.

4. STORAGE

Keep out of the sight and reach of children. Do not store above 25°C. Store in original package. Do not use ORALSAN OxyGel after the expiry date stated on the label or carton.

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