ORALSAN®

Wf**3**

Compact
Water Flosser

Pag 2

Pag 9

Pag 16

www.idsdental.it


Oralsan® wf3 - Compact Water Flosser GB

Thanks for your purchase of Oralsan® water flosser wf3. Please read this user manual carefully and keep it in a safe place for future reference. If you need any help, please contact us.

Model.....	WF3
Input	5V 1A
Material	ABS PC
Battery capacity.....	2500 mAh
Potenza di esercizio	21W
Water pressure range	40-110psi
Charging time	3-5 hours
Waterproof level	IPX7
Reservoir capacity	200 ml
Water flow method	Pulsating water

Product parts diagram

1. Nozzle tip
2. Power button
3. Mode selection button
4. Mode indicator lights
5. Charging indicator light
6. Nozzle release button
7. Charging port
8. Water tank
9. Lid for water tank

Accessories

USB charging
cable

Tongue
cleaner tip

Orthodontic
tip

Periodontal
tip

Jet
tip

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING.

DANGER:

To reduce the risk of electrocution:

Do not handle charger with wet hands.

Do not place charger in or drop charger into water or other liquid.

Check the charger cord for damage before the first use and during the life of the product.

Use this product only as indicated in these instructions or as recommended by your dental professional.

Only use tips and accessories that are recommended by H2ofloss.

Do not drop or insert any foreign object into any opening on the product.

Keep charging unit and cord away from heated surfaces.

Do not use this product outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.

Do not use iodine, saline (salt solution), or water insoluble concentrated essential oils in this product. Use of these can reduce performance and will shorten the life of the product. Use of these products may void warranty.

Fill reservoir with water or other dental professional recommended solution only.

Remove any oral jewelry prior to use of this product.

Do not use if you have an open wound on your tongue or in your mouth.

Maintenance

Make sure to throw out the water inside the tank, and run the appliance for 2 to 3 seconds to drain the water after usage. Otherwise bacteria could collect in the remaining water. Do not submerge the unit in water or any other liquid. Wipe the exterior when necessary with a soft cloth and mild non-abrasive cleanser. Reservoir may be washed with soapy water.

Wash the dirt off from inside using a brush or similar, and wipe off with a cloth. Dry the reservoir when not using for a long period.

Battery Care and Disposal

The unit contains a Lithium battery.

At the end of the battery life, please recycle the unit at your nearest recycle center per local or state requirements. The battery is not replaceable.

Accessories

5 x Tips (see box)

1 x User manual

1 x USB Power Cable

1 x Travel case

Warranty

WF3 will repair or replace (at our option) your unit free of charge for 2 months from the date of purchase if the appliance is defective in workmanship or materials. This guarantee does not cover: damage due to improper use, normal wear or use as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorized persons and if original WF3 parts are not used. To obtain service within the guarantee period, hand in or send the complete appliance with your purchase receipt to an authorized WF3 distributor/customer service center in your country.

ELCAT S.N.C - VIA P. CRESPI, 10 - MILANO

Tel. +39 02 26111520 - fax +39 02 26143958 - e-mail: info@elcatmilano.com

Purchaser is responsible for all postage and freight charges and any other expenses incidental to claiming a gainst this Warranty.

Note: Excessive heat will cause machine to stop when in use, Switch off machine, cool for 2 hours and re-start

Made in China
Printed in China

WARNING:

To reduce the risk of burns, electrocution, fire, or injury to persons:

Do not plug this device into a voltage system that is different from the voltage system specified on the device or charger.

improper voltage will damage product and VOID warranty.

If product is opened/disassembled for ANY reason, warranty is VOID.

Do not use this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.

Do not direct water under the tongue, into the ear, nose or other delicate areas. This product is capable of producing pressures that may cause serious damage in these areas.

Do not direct water into the nose or ear. The potentially deadly amoeba, *Naegleria fowleri*, may be present in some tap water or unchlorinated well water and may be fatal if directed into these areas.


If your physician advised you to receive antibiotic premedication before dental procedures, you should consult your dentist before using this instrument or any other oral hygiene aid.

Close supervision is necessary when this product is used by, on or near children or invalids.

Children should be supervised to ensure that they do not play with the appliance.

SAVE THESE INSTRUCTIONS.


Charging Unit Prior to First Use

The unit is charged through devices with 5V USB's power capabilities such as mobile phone, laptop, USB power wall charger, etc. The charge indicator on the handle will flash when charging and will stay on continuously when the unit is fully charged. A complete recharge will take 3-5 hours.

Filling the Reservoir


Lift flip top on reservoir and fill the reservoir with lukewarm water.

Inserting and Removing Tips

Insert tip into the center of the knob at the top of the Water Flosser handle. Press firmly until the tip clicks into place. To remove the tip from the handle, press the tip eject button and pull tip from the handle.

CAUTION: Do not press tip eject button while unit is running.

Adjusting the Pressure Setting

You can adjust the pressure setting on your Water Flosser by pressing the mode button on the handle. Start with the lowest setting and increase to your liking.

Pressure
Setting

When Finished

Turn the unit off. Empty any liquid left in the reservoir.

TIP USAGE

Jet tip

Been over the sink and guide the jet to your teeth and gums. Glide the handle switch up to start the water stream. Keep your mouth slightly open to allow the water to run into the sink. Direct the stream of water towards the gumline at 90 degrees, so that contact is made where the tooth and gums meet.


Manoeuvre the flow of water from tooth to tooth remembering to pause briefly between each tooth.

Periodontal tip

Place the periodontal tip in contact between teeth and the gumline. Lean the soft tip against a tooth at a 45-degree angle and gently place the tip under the gumline, into the pocket. Start the irrigator by setting it at the lowest pressure and clean by gently following gumline, inserting the periodontal tip in between teeth. Continue tracing along the gumline to massage the gums and safely deliver solution into the gingival pockets below the gumline.

Periodontal tip

Place the periodontal tip in contact between teeth and the gumline. Lean the soft tip against a tooth at a 45-degree angle and gently place the tip under the gumline, into the pocket. Start the irrigator by setting it at the lowest pressure and clean by gently following gumline, inserting the periodontal tip in between teeth. Continue tracing along the gumline to massage the gums and safely deliver solution into the gingival pockets below the gumline.

Tongue cleaner tip

Lean your head over the sink. Glide the pressure button to the lowest setting. Start the irrigator and place the tongue cleaner at the back of your tongue and move it forward to clean your tongue. Repeat to cover entire tongue. Set the pressure to the lowest setting and adjust to a higher setting for more water pressure as it feels comfortable.

Troubleshooting

Problem	Cause	Corrective action
Reservoir leaks	Seal ring is abnormal	Remove the seal ring and reinstall
Low water pressure	Low battery	Recharge the device
	Nozzle head blocks	Replace with new nozzle heads
Lasts only few minutes with full charge	Dead battery	Contact retailers
	Charging cord	Try with new charging cable
Motor makes sound but no water coming out	Nozzle head blocks	Replace with new tips
	Suction valve dislocation	Keep device on, fill water to the hole where nozzle head goes into
Motor doesn't make sound	Low Battery (Keep the device idle for more than 3 month)	Activate the battery by recharging it for 16 hours
Product stopped working suddenly	Overheating protection is activated (Continuous use 5 tanks of water with highest pressure)	Restart the device in 2 hour

If your problem still hasn't been fixed by above troubleshooting, please kindly contact with the retailers for help.

- It may occur below situation when first time use this machine or left it unused for a long time: The nozzle head won't spray water immediately after turning on the machine. In this case, please follow below instructions to make it work:

Add proper amount of water into the place where tips are inserted in.

Insert the tip in and press power button to turn on the device.

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